

Surgical special protect apron

Designed for wearing under the gown.

Releasing the velcro-fastener allows the apron to drop without removing the sterile gown.



Lead material
Light: light lead rubber material

Lead equivalent
front 0,50 mm Pb / back 0,25 mm Pb
(upon request: front 0,25; 0,35 mm Pb /
back 0,25 mm Pb)

Standard closure
Velcro

Colours
Standard: royalblue

Size

Medium M : Hip 60 / Shoulder 48

Lenght cm:

90 130090

100 130100

110 130110

120 130120

130 130130

Large L: Hip 75 / Shoulder 53

Lenght cm:

100 140100

110 140110

120 140120

130 140130

Alterations reserved !

EWEM-JANUS · Inh. Waltraud Janus · X-ray accessories

D-94405 Landau/Isar · Zum Alten Hofgarten 22 · phone +49 (0) 9951-8051 · fax +49 (0) 9951-6399

Internet: <http://www.ewem-janus.de> · e-mail: info@ewem-janus.de